

Beetle Lysis-Juice (1 02 512)

Components included:

Beetle Lysis-Juice	100ml Buffer for Screening with Beetle-Luciferase – without substrate. Store at +4°C.
D-Luciferin	1 vial (A)→ for 100 ml Beetle Lysis-Juice p. a. free acid, firefly (synth.). Store at -20°C.
ATP	1 vial (B)→for 100 ml Beetle Lysis-Juice Adenosintriphosphat, Store at +4°C.

Reconstruction:

Beside the black buffer-bottle the kit is equipped with two further tubes which are labelled with **A** D-Luciferine and **B** ATP. D-Luciferine and ATP are white powders which have to be dissolved and mixed into the buffer by shaking gently. Afterwards the finished measuring buffer can be aliquoted and long-term stored at -80°C.


Note: Do not combine this buffer with marine luciferases because the buffer contains triton. The combination of triton with the substrate of marine luciferases (Coelenterazine) leads to autoluminescence.

Standard Protocol

Program Luminometer:


For the measuring we suggest measurement duration of 5 sec. and a delay of 2 sec. for each sample.

- 1.) Add 100 µl of **Beetle Lysis-Juice** to the cells (free of medium).
- 2.) To make sure the lysis was completed please mix the dissolved cells 3 times in the pipette tip.
- 3.) After that incubate the solution at least for 5 min. (up to 20 min. possible).
- 4.) Measurement starts after 2 sec. delay for 5 sec. duration.



Beetle Lysis-Juice (1 02 512)

Standard procedure:



Optional: parallel measurement in Duo- or Triple-Well System (Firefly/ Renilla/ Gaussia/AP/ β Gal)

Split the Cell Lysate and transfer to different tubes for measuring several luciferases in one sample.

